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Dear Valued Streck Customer,

This letter is to notify you of updates that have been made to the Instructions for Use (IFU) document for Cell-Free DNA BCT<sup>®</sup> CE.

The Plasma Isolation section has been updated with a single spin protocol that matches the protocol on our U.S. IVD cleared assay. This was done to reduce customer confusion about which spin protocol to use. The following note was also added to this section: "Note: Cell-Free DNA BCT has been validated for a maximum centrifugation of 3,000 x g for 10 minutes. Exceeding these limits may result in breakage."

The References section was updated with the following statement: "2. ISO 6710, Single-use containers for human venous blood specimen collection." This statement was added so that customers know which standard the glass tubes have been validated to.

Cell-Free DNA BCT CE is for Export Only. Not for use in the U.S.

Streck is committed to assisting our customers with this document change. If you have any questions, please contact Streck Technical Services at 402-691-7510 or 800-843-0912, or a Sales Representative at 402-333-1982 or 800-843-0912.

Kind Regards,

A handwritten signature in black ink that reads "Bill Rock". The signature is fluid and cursive, with the first letters of "Bill" and "Rock" being capitalized and prominent.

Bill Rock  
Product Manager  
[brock@streck.com](mailto:brock@streck.com)

## INSTRUCTIONS FOR USE

Cell-Free DNA BCT® CE is a direct draw whole blood collection tube intended for collection, transport and storage of blood samples. **This product is FOR EXPORT ONLY, not to be sold in the United States.**

## SUMMARY AND PRINCIPLES

Cell-Free DNA BCT CE stabilizes cell-free plasma DNA as well as preserves cellular genomic DNA present in nucleated blood cells and circulating epithelial cells (tumor cells) found in whole blood.

Accurate analysis of cf-DNA can be compromised by sample handling, shipping and processing, causing lysis of nucleated blood cells and subsequent release of cellular genomic DNA. Additionally, degradation of cf-DNA due to nuclease activity can be problematic.

The preservative reagent contained in Cell-Free DNA BCT CE stabilizes nucleated blood cells, preventing the release of cellular genomic DNA, and inhibits nuclease mediated degradation of cf-DNA, contributing to the overall stabilization of cf-DNA. Samples collected in Cell-Free DNA BCT CE are stable for up to 14 days at temperatures between 6 °C to 37 °C, allowing convenient sample collection, transport and storage.

The preservative reagent contained in Cell-Free DNA BCT CE stabilizes circulating epithelial cells (tumor cells) in whole blood for up to 7 days at temperatures between 15 °C to 30 °C. Cell-Free DNA BCT CE is for use by laboratory professionals.

## REAGENTS

Cell-Free DNA BCT CE contains the anticoagulant K<sub>2</sub>EDTA and a cell preservative in a liquid medium. No additional materials are provided or required.

## PRECAUTIONS

- For In Vitro Diagnostic Use.
- Do not freeze specimens collected in glass Cell-Free DNA BCT CE.
- Do not use tubes after expiration date.
- Do not use tubes for collection of materials to be injected into patients.
- Product is intended for use as supplied. Do not dilute or add other components to Cell-Free DNA BCT CE.
- Overfilling or underfilling of tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.

### CAUTION

- Glass has the potential for breakage; precautionary measures should be taken during handling of glass tubes.
  - All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection. Dispose of in accordance with federal, state and local regulations. Avoid contact with skin and mucous membranes.
  - Product should be disposed with infectious medical waste.
  - Remove stopper by either gently rocking the stopper from side to side or by grasping with a simultaneous twisting and pulling action. A "thumb roll" procedure for stopper removal is not recommended, as tube breakage and injury may result. Reinsert stopper by gently pushing stopper onto tube with a simultaneous twisting action.
7. SDS can be obtained at [streck.com](http://streck.com) or by calling 1-402-691-7510.

## STORAGE AND STABILITY

- When stored at 2 °C to 30 °C, empty Cell-Free DNA BCT CE is stable through expiration date.
- Short-term storage at 2 °C to 40 °C is acceptable for empty Cell-Free DNA BCT CE for up to 14 days.
- Do not freeze empty Cell-Free DNA BCT CE. Proper insulation may be required for shipment during extreme temperature conditions.
- Sample storage/stability:

	Sample Type		
	Cell-Free DNA	Cellular Genomic DNA	Epithelial Cells (Tumor Cells)
Sample Stability	14 days	14 days	7 days
Sample Storage Temperature	6 °C to 37 °C	6 °C to 37 °C	15 °C to 30 °C


## INDICATIONS OF PRODUCT DETERIORATION

- Cloudiness or precipitate visible in reagent of empty tube.
- If indications of product deterioration occur, contact Streck Technical Services at 1-402-691-7510 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

## INSTRUCTIONS FOR USE

For a video demonstration, visit [streck.com/mixing](http://streck.com/mixing).

- Collect specimen by venipuncture according to JCLLS GP4-A3'.  
**Prevention of Backflow** - Since Cell-Free DNA BCT CE contains chemical additives, it is important to avoid possible backflow from the tube.  
To guard against backflow, observe the following precautions:
  - Keep patient's arm in the downward position during the collection procedure.
  - Hold the tube with the stopper in the uppermost position so that the tube contents do not touch the stopper or the end of the needle during sample collection.
- Cell-Free DNA BCT CE should be drawn after the EDTA tube and before the fluoride oxalate (glycolytic inhibitor) tube. If a Cell-Free DNA BCT CE tube immediately follows a heparin tube in the draw order, Streck recommends collecting a non-additive or EDTA tube as a waste tube prior to collection in the Cell-Free DNA BCT CE.
- Fill tube completely.
- Remove tube from adapter and immediately mix by gentle inversion 8 to 10 times. Inadequate or delayed mixing may result in incorrect analytical results or poor product performance. One inversion is a complete turn of the wrist, 180 degrees, and back per the figure below:

- After collection, transport and store tubes within the recommended temperature range.

Report any serious incident to Streck and appropriate regulatory entities including the competent authority of the member state in which the user and/or patient is established as applicable.

## Note:

- For best results, a 21G or 22G needle is advised. Slower fill times may be observed when using a smaller gauge needle.
- When using a winged (butterfly) collection set for venipuncture and the Streck Cell-Free DNA BCT CE is the first tube drawn, a non-additive or EDTA discard tube should be partially drawn first in order to eliminate air or "dead space" from the tubing.
- Cell-Free DNA BCT CE does not dilute blood samples; therefore, no dilution factor correction is necessary.
- As in the case with most clinical laboratory specimens, hemolysis, icterus and lipemia may affect the results obtained on blood samples preserved with Cell-Free DNA BCT CE.

## PLASMA ISOLATION

- Centrifuge blood collection tubes containing whole blood for 10 minutes at room temperature at 1,600 x g. Carefully unload the tubes from the centrifuge without disturbing the buffy coat.
- Pipette maximum 5ml plasma from the specimen tube to a labeled 15ml centrifuge tube.
- Centrifuge plasma in 15ml centrifuge tube for 10 minutes at room temperature at 3,220 x g. Note: Centrifuge speed may be adjusted down to 3,200 x g, if centrifuge does not support 3,220 x g.
- Pour contents of 15ml centrifuge tube for each specimen into an appropriately labeled 5ml conical screw cap tube.
- Isolate cell-free DNA using extraction kit manufacturer's instructions including any adaptations as required.

Note: Cell-Free DNA BCT has been validated for a maximum centrifugation of 3,000 x g for 10 minutes. Exceeding these limits may result in breakage.

For optimal results for all of the above protocols, include a Proteinase K treatment step ( $\geq 30$  mAU/mL digest) at 60 °C in the presence of chaotropic salts for 1 hour when extracting cell-free DNA.

## Cellular Genomic DNA

- To separate the white blood cells, either lyse the red blood cells and wash, or centrifuge whole blood and collect the buffy coat layer.
- Isolate genomic DNA per kit manufacturer instructions.

For optimal results, include a Proteinase K treatment step ( $\geq 30$  mAU/ml digest) at 60 °C in the presence of chaotropic salts for 2 hours when extracting cellular genomic DNA.

## FREEZING AND THAWING

### PLASMA

- To Freeze: For long-term storage, after spinning, collect and transfer the upper plasma layer to a cryogenic tube (not provided) and freeze at -20 °C or -80 °C.
- To Thaw: Thaw cryogenic tubes at appropriate temperature as specified in your protocol.  
**Note:** If cryoprecipitates form in the plasma, vortex the tube for 30 seconds after thawing. Do not centrifuge the plasma.

## LIMITATIONS

- For single use only.
- Samples drawn in other anticoagulants or preservatives may cause coagulation in Cell-Free DNA BCT CE.
- Specimen transport via pneumatic tube system is not advised.

## REFERENCES

- JCLLS Standard Phlebotomy Guideline GP4-A3.
- ISO 6710, Single-use containers for human venous blood specimen collection

## ORDERING INFORMATION

Please call our Customer Service Department at 1-402-691-7510 for assistance. Additional information can be found online at [streck.com](http://streck.com).

## GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [streck.com](http://streck.com).

## CHANGES FROM PREVIOUS VERSION

Updated Plasma Isolation protocol and References.

Australia Patent AU2003254755  
 Canada Patent CA2,917,912  
 Europe Patent EP2228453B1; EP2626438A1; EP2814981; EP1816461  
 Germany Patent DE60201322817.5; DE202010048559  
 United States Patent US9,657,227; US9,926,590; US10,144,955; US10,294,513; US10,091,984  
 Others Pending  
 See [streck.com/patents](http://streck.com/patents) for patents that may be applicable to this product.

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