

[コストパフォーマンスに優れた次世代磁気ビーズをChIP-seqで使用 研究者の声【43】](#)

東京大学 山中 総一郎 先生

プロトコル紹介

>> [ChIP-seq 実験にベリタス取扱製品を使用 研究者の声【16】](#) にて紹介した既存プロトコルを [DynaGreen版](#) に一部編集 (2024/09/18)

【実際に行ったDynaGreenを利用した ChIP-seq のプロトコル】

1. 細胞の固定

↓マウス胚性幹細胞を **40,000** 細胞分取する (in low-binding tube)

↓remove supernatant

↓suspend with 1 ml DMEM (10% FBS)

↓add 67 μ l, 16% Formaldehyde (Pierce) (final 1%)

↓RT, 10min, rotation

↓add 62.5 μ l 2 M glycine (final 125 mM)

↓RT, 5 min, rotation

↓centrifuge, 6,000 rpm, 1 min x2

↓take supernatant (~900 μ l)

↓centrifuge, 6,000 rpm, 1 min

↓take supernatant completely > stock in -80°C freezer

2. DynaGreen のブロッッキング

↓take **25 μ l DynaGreen beads** into two tubes per sample (one for IP and one for pre-clear)

↓wash with 800 μ l, 0.5% BSA in PBS, twice

↓suspend with 1 ml, 0.5% BSA in PBS

↓4°C, rotation, 1 h

3. IP 用 DynaGreen の調製

↓add 2 μ g anti-body (H3K9me3)

↓4°C, rotation, o/n

4. ソニケーション

Frozen cells

↓add 1 ml Swelling Buffer, 4°C, rotation, 20 min

↓centrifuge, 8,000 rpm, 1 min x2

↓take supernatant (~900 μ l)

↓centrifuge, 8,000 rpm, 1 min

↓take supernatant completely

↓suspend with 100 μ l 1x Sonication buffer (Covaris)

↓sonication [Branson: 30%, total ON time=2 min, ON/OFF=2 sec/3 sec]

RPDBM-24-1428

(2024.9.18)

- ↓centrifuge, 13,500 rpm, 5 min
- ↓save supernatant to low-binding tube
- ↓centrifuge, 13,500 rpm, 10 min
- ↓save supernatant to low-binding tube
- ↓dilute with 1 ml of Abcam RIPA buffer

5. pre-clear

2 で調製した“pre-clear beads”

- ↓wash with 1ml, Abcam RIPA buffer x2 回
- ↓suspend with 20 µl Abcam RIPA buffer
- ↓add to cell extract
- ↓4°C, rotation, 1 h
- ↓save supernatant
- ↓aliquot 100 µl for [WCE] > the rest of the cell extract(~1 ml) is for IP

6. IP

3 で調製した“beads conjugated with anti-body”

- ↓wash with 1 ml, Abcam RIPA buffer x4 回
- ↓suspend with 20 µl Abcam RIPA buffer
- ↓add beads to 1 ml cell extract
- ↓4°C, rotation, o/n

7. Beads wash > elution

IP を行った beads

- ↓wash with 1 ml, LOW buffer x3 回
- ↓wash with 1 ml, HIGH buffer x3 回
- ↓add 200 µl direct elution buffer
- ↓incubate at 65°C, 15 min, with constant vortexing
- ↓save supernatant

WCE_100 µl of WCE

- ↓add 100 µl direct elution buffer
- ↓incubate at 65°C, 15 min, with constant vortexing

8. ProK treatment > reverse crosslinking

IPed sample and WCE sample

- ↓transfer to screw tube
- ↓add 5 µl proteinase K (Roche)
- ↓incubate 37 °C, >6 h
- ↓65°C, >6 h

RPDBM-24-1428

(2024.9.18)

9. DNA 回収

↓add 20 µl 3 M NaOAc, 200 µl Phenol:Chloroform (Thermo)

↓centrifuge, 13,500 rpm, 10 min

↓save supernatant to low-binding tube

↓+1 µl pellet paint (Merck)

↓-80°C, 15 min

↓centrifuge, 14,000 rpm, 30 min

↓remove supernatant

↓add 150 µl, 70% ethanol

↓centrifuge, 14,000rpm, 5 min

↓remove supernatant

↓dry up RT, 10 min

suspend with 35 µl 0.1M Tris-HCl.

10. ChIP-seq library preparation

Check the quality of IPed sample with qPCR

↓prepare the DNA library for sequencing with QIAGEN QIAseq Ultralow Input library Kit

According to the manufacturer's instruction

↓check the DNA size with Bioanalyzer

使用した Buffer の組成

- Swelling buffer

20 mM HEPES(pH7.9)

1.5 mM MgCl₂

10 mM KCl

0.1% NP-40

1mM DTT (add before use)

- Abcam RIPA buffer

(RIPA buffer in Abcam protocol)

50 mM Tris-HCl, pH8.0

150 mM NaCl

2 mM EDTA pH8

1% NP-40

0.5% Sodium Deoxycholate (DOC)

0.1% SDS

RPDBM-24-1428

(2024.9.18)

- Low buffer

(Wash buffer in Abcam protocol)

0.1% SDS

1% Triton X-100

2 mM EDTA pH8

150 mM NaCl

20 mM Tris-HCl, pH8.0

- High buffer

(Final wash buffer in Abcam protocol)

0.1% SDS

1% Triton X-100

2 mM EDTA pH8

500 mM NaCl

20 mM Tris-HCl, pH8.0

- Direct elution buffer

10 mM Tris-HCl, pH8.0

5 mM EDTA

300 mM NaCl

0.5% SDS

株式会社ベリタス 〒105-0013 東京都港区浜松町1-18-16 住友浜松町ビル6階

TEL 03-5776-0078 FAX 03-5776-0076

技術的なお問い合わせは：TEL 03-5776-0040 E-mail Tech_support@veritastk.co.jp

RPDBM-24-1428

(2024.9.18)